

Baby Head Band™ and Circuit Bumpers™

A more comfortable and versatile way to
secure nasal prongs

respiralogics™

Baby Head Band and Circuit Bumpers are intended to secure the nasal interface during oxygen therapy, nCPAP or NIV therapy.

The Baby Head Band is made of breathable, non-latex foam that absorbs and wicks away moisture leaving the scalp and skin dry. The soft material of the head band allows it to conform to even the tiniest of heads while maintaining the circuit position. The soft outer covering of the Baby Head Band allows for head movement while holding the nasal prongs and circuit stable.

The new reversible Circuit Bumpers secure and stabilize the position of the inspiratory and expiratory limbs to the soft outer side of Head Band while providing optimal prong positioning for infants of all shapes and sizes.


The Baby Head Band and Circuit Bumpers work in conjunction with neonatal/pediatric ventilator circuits and oxygen therapy tubing. Designed for use in NICU, PICU, surgery and clinics for stable positioning of nasal interface.

Ordering Information

30501	Baby Head Band and Circuit Bumpers – Size Extra Small (18 to 23 cm)	15 / box
30502	Baby Head Band and Circuit Bumpers – Size Small (24 to 30 cm)	15 / box
30503	Baby Head Band and Circuit Bumpers – Size Medium (31 to 36 cm)	15 / box
30504	Baby Head Band and Circuit Bumpers – Size Large (36 to 42 cm)	15 / box

Companion Products

1002	Babi.Plus® nCPAP Kit	B020	Babi.Plus nCPAP Nasal Cannula, Sz 0
1010	Babi.Plus nCPAP Kit, Small	B021	Babi.Plus nCPAP Nasal Cannula, Sz 1
1020	Babi.Plus nCPAP Kit, Large	B022	Babi.Plus nCPAP Nasal Cannula, Sz 2
1005	Venti.Plus™ Single Pole Clamp	B023	Babi.Plus nCPAP Nasal Cannula, Sz 3
1006	Venti.Plus Dual Pole Clamp	B024	Babi.Plus nCPAP Nasal Cannula, Sz 4
30100	Baby Nose Bumper™ and RespiraGel Mustache™	B025	Babi.Plus nCPAP Nasal Cannula, Sz 5
30150	Baby Nose Bumper Replacement	B026	Babi.Plus nCPAP Nasal Cannula, Sz 6
30400	Baby Chin Strap™	B027	Babi.Plus nCPAP Nasal Cannula, Sz 7
30510	Baby Band™ for Chin and Head		



FEATURES AND BENEFITS

Fashioned to provide a comfortable fit for all infants

Sized by head circumference for easy and quick application

Reversible hook-and-loop bumpers allow secure attachment and minimize pressure points

Allows caregivers to freely move infant without interrupting therapy

Opened at the crown for fontanel checks and to accommodate scalp vein IVs

Available in four sizes for very low birth weight up to 10 kg infants

Single patient use, non-latex material

Respiralogics products are available from Respiralogics and our specialty respiratory care distributors. For our latest product information visit www.respiralogics.com or email us at 4info@respiralogics.com.



Respiralogics is a Global Respiratory Solutions, Inc. Company
3545 Airway Drive, Suite 104 • Reno, NV 89511 • 775 954 0160 • www.respiralogics.com

© 2019 Respiralogics. All rights reserved. Respiralogics, the Respiralogics logo, Baby Head Band and Circuit Bumpers, Baby Nose Bumper and RespiraGel Mustache, and Baby Band for Chin and Head are trademarks of Global Respiratory Solutions, Inc. Babi.Plus is a registered trademark and Venti.Plus is a trademark of GaleMed Inc. C30500.2