

Babi.Plus® Bubble PAP Valve

respiralogics™

Accurate system for delivery of BCPAP therapy

The Babi.Plus Bubble PAP Valve provides a convenient method of delivering continuous PAP (positive airway pressure) to infants weighing less than 10 kg. The patent-pending design delivers accuracy and stability throughout the infant's course of therapy, freeing the clinician to focus on patient care.

The Babi.Plus Bubble PAP Valve was the first to receive FDA 510k market clearance for delivery of positive end expiratory pressure in patients weighing < 10 Kg in hospital environments to increase end lung pressure above atmospheric in a constant flow condition. The PAP is adjustable from 1 to 10 cm H₂O, using gas flows from 1 to 12 liters per minute for precise pressure delivery. The Bubble PAP Valve's ergonomic design allows pressure to be easily set without the cumbersome and time consuming tasks normally associated with bubble devices.

The Babi.Plus Bubble PAP Valve jar has an integrated ISO compliant equipment mount, which makes it easy to secure the jar onto a standard equipment pole with a clamp. The Venti.Plus™ Single or Dual Pole Clamp companion products are designed to pair with Babi.Plus products and also work well with other products requiring pole clamps.

The Babi.Plus Bubble PAP Valve can be used in all neonatal critical care units, delivery rooms, pediatric critical care units and special procedure units.

Ordering Information

1002	Babi.Plus BCPAP Kit, includes both BCPAP Valve and Pressure Regulator	12/case
1000	Babi.Plus Bubble PAP Valve 0-10 cm H ₂ O	12/case

Companion Products

1010	Babi.Plus nCPAP Kit, Small	30100	Baby Nose Bumper™ and RespiraGel Mustache™
1020	Babi.Plus nCPAP Kit, Large	30150	Baby Nose Bumper Replacement
1030	Babi.Plus Pressure Relief Manifold	30301	Baby Cap™ & Circuit Bumpers™ Size 1
B020	Babi.Plus nCPAP Nasal Cannula, Sz 0	30302	Baby Cap & Circuit Bumpers, Size 2
B021	Babi.Plus nCPAP Nasal Cannula, Sz 1	30303	Baby Cap & Circuit Bumpers, Size 3
B022	Babi.Plus nCPAP Nasal Cannula, Sz 2	30304	Baby Cap & Circuit Bumpers, Size 4
B023	Babi.Plus nCPAP Nasal Cannula, Sz 3	30305	Baby Cap & Circuit Bumpers, Size 5
B024	Babi.Plus nCPAP Nasal Cannula, Sz 4	30306	Baby Cap & Circuit Bumpers, Size 6
B025	Babi.Plus nCPAP Nasal Cannula, Sz 5	30307	Baby Cap & Circuit Bumpers, Size 7
B026	Babi.Plus nCPAP Nasal Cannula, Sz 6	30400	Baby Chin Strap™
B027	Babi.Plus nCPAP Nasal Cannula, Sz 7	30501	Baby Head Band™, Size Extra Small
GB10	GiO™ 1 Pressure Gauge	30502	Baby Head Band, Size Small
1005	Venti.Plus Single Pole Clamp	30503	Baby Head Band, Size Medium
1006	Venti.Plus Dual Pole Clamp	30504	Baby Head Band, Size Large
		30510	Baby Band™ for Chin and Head

Respiralogics products are available from Respiralogics and our specialty respiratory care distributors. For our latest product information visit www.respiralogics.com or email us at 4info@respiralogics.com.



Respiralogics is a Global Respiratory Solutions, Inc. Company
3545 Airway Drive, Suite 104 • Reno, NV 89511 • 775 954 0160 • www.respiralogics.com

© 2019 Respiralogics. All rights reserved. Respiralogics and the Respiralogics logo, Baby Cap and Circuit Bumpers, Baby Nose Bumper and RespiraGel Mustache, Baby Chin Strap, Baby Head Band and Circuit Bumpers, and Baby Bands for Chin and Head are trademarks of Global Respiratory Solutions, Inc. Babi.Plus is a registered trademark and GiO and Venti.Plus are trademarks of Galemel Corporation. C1000.5



FEATURES AND BENEFITS

Accurate method to deliver PAP for infants weighing ≤10 kg

Rotating pressure lever is adjustable throughout clinical flow rate and pressure ranges

Pressure settings clearly marked for precise pressure adjustment

Maximum and Minimum fluid lines clearly visible

Positive engagement pressure lever reduces risk of inadvertent changes

360° swivel inlet connector accepts 15 mm OD by 22 mm ID expiratory circuit connector

Water level adjustment port allows fluid to be added or removed without disconnect

Gas discharge ports direct expiratory gas away from the clinician

Latex free, hypoallergenic material

Single-patient use closed system minimizes the incidence of atmospheric fluid contamination

Ready for immediate set up in all clinical environments