

Venti.Plus™ Single Pole Clamp and Dual Pole Clamp

respiralogics™

Convenient method to attach Babi.Plus Bubble PAP Valve

The Venti.Plus Single Pole Clamp and Venti.Plus Dual Pole Clamp are sturdy pole clamps designed to support the Babi.Plus® Bubble PAP Valve 0–10 cm H₂O on a typical equipment pole. The Venti.Plus Single Pole Clamp is durable and intended to support weight of ≤ 3.0 kg. The unique Venti.Plus Dual Pole Clamp is intended to support a tandem weight of ≤10 kg of two products on an equipment pole, specifically to hold the Bubble PAP Valve and a humidifier. The Pole Clamps provide a convenient method to mount the Babi.Plus Bubble PAP Valve to a variety of equipment poles.

The Venti.Plus Pole Clamps are ISO compliant for mounting on poles 2.0 to 3.2 cm in diameter. The Venti.Plus Pole Clamps are a companion product for the Babi.Plus Bubble PAP Valve and allow for mounting of the Bubble PAP Valve, humidifiers or other devices.


Ordering Information

1005	Venti.Plus Single Pole Clamp	1 each
1006	Venti.Plus Dual Pole Clamp	1 each

Companion Products

1002	Babi.Plus Bubble CPAP Kit
1000	Babi.Plus Bubble PAP Valve 0-10 cm H ₂ O
1016	Babi.Plus nCPAP Complete Kit, Small
1026	Babi.Plus nCPAP Complete Kit, Large

Respiralogics products are available from Respiralogics and our specialty respiratory care distributors. For our latest product information visit www.respiralogics.com or email us at 4info@respiralogics.com.



FEATURES AND BENEFITS

- Specifically designed to hold and use with the Babi.Plus Bubble PAP Valve and humidifiers
- Easy to use for attachment to typical equipment pole
- Easy application to equipment pole with thumb screw mechanism
- Sturdy clamp that may be used with multiple devices
- May be cleaned easily between clinical uses

Respiralogics is a Global Respiratory Solutions, Inc. Company
3545 Airway Drive, Suite 104 • Reno, NV 89511 • 775 954 0160 • www.respiralogics.com

© 2019 Respiralogics. All rights reserved. Respiralogics and the Respiralogics logo are trademarks of Global Respiratory Solutions, Inc. Babi.Plus is a registered trademark and Venti.Plus is a trademark of Galemed Corporation. C1005.5