

# Babi.Plus® Pressure Relief Manifold

Prevent pressure build up in the circuit

respiralogics™

12.5 cm H<sub>2</sub>O Pressure Relief Manifold with Gas Supply Tube and Oxygen Sensor Port is a single use device for use with patients weighing ≤ 10 kg. The manifold is designed to prevent breathing circuit pressure from exceeding the preset pressure limit.

The Babi.Plus Pressure Relief Manifold also includes a Gas Supply Tube, an Oxygen Sensor Port and capped port for an oxygen sensor analyzer.

The Babi.Plus Pressure Relief Manifold is designed to pair with the Babi.Plus Bubble PAP Valve 0-10 cm H<sub>2</sub>O and Babi.Plus Cannula Kits companion products. The Babi.Plus Pressure Relief Manifold can be used in neonatal critical care units, delivery rooms, pediatric critical care units and special procedure units.

## Ordering Information

1030	Babi.Plus Pressure Relief Manifold	12/case
1002	Babi.Plus BCPAP Kit	12/case

## Companion Products

1000	Babi.Plus Bubble PAP Valve	30100	Baby Nose Bumper™ and RespiraGel Mustache™
1010	Babi.Plus nCPAP Kit, Small	30150	Baby Nose Bumper Replacement
1020	Babi.Plus nCPAP Kit, Large	30301	Baby Cap & Circuit Bumpers™, Size 1
1030	Babi.Plus Pressure Relief Manifold	30302	Baby Cap & Circuit Bumpers, Size 2
B020	Babi.Plus nCPAP Nasal Cannula, Sz 0	30303	Baby Cap & Circuit Bumpers, Size 3
B021	Babi.Plus nCPAP Nasal Cannula, Sz 1	30304	Baby Cap & Circuit Bumpers, Size 4
B022	Babi.Plus nCPAP Nasal Cannula, Sz 2	30305	Baby Cap & Circuit Bumpers, Size 5
B023	Babi.Plus nCPAP Nasal Cannula, Sz 3	30306	Baby Cap & Circuit Bumpers, Size 6
B024	Babi.Plus nCPAP Nasal Cannula, Sz 4	30307	Baby Cap & Circuit Bumpers, Size 7
B025	Babi.Plus nCPAP Nasal Cannula, Sz 5	30400	Baby Chin Strap™
B026	Babi.Plus nCPAP Nasal Cannula, Sz 6	30501	Baby Head Band™, Size Extra Small
B027	Babi.Plus nCPAP Nasal Cannula, Sz 7	30502	Baby Head Band, Size Small
GB10	GiO™ 1 Pressure Gauge	30503	Baby Head Band, Size Medium
1005	Venti.Plus™ Single Pole Clamp	30504	Baby Head Band, Size Large
1006	Venti.Plus Dual Pole Clamp	30510	Baby Band™ for Chin and Head

Respiralogics products are available from Respiralogics and our specialty respiratory care distributors. For our latest product information visit [www.respiralogics.com](http://www.respiralogics.com) or email us at [4info@respiralogics.com](mailto:4info@respiralogics.com).




**FEATURES AND BENEFITS**

- Easy to use manifold for delivery of non-invasive therapy for infants weighing ≤ 10 kg
- 15 mm ID x 22 mm OD inlet adaptor interfaces with all humidifier chambers
- Integrated oxygen supply tubing with fits all connector for mixed gas delivery
- Constructed of Phthalate Free, Latex free and hypoallergenic material
- Pre-assembled, single patient use manifold

Respiralogics is a Global Respiratory Solutions, Inc. Company  
3545 Airway Drive, Suite 104 • Reno, NV 89511 • 775 954 0160 • [www.respiralogics.com](http://www.respiralogics.com)

© 2019 Respiralogics. All rights reserved. Respiralogics and the Respiralogics logo, Baby Cap and Circuit Bumpers, Baby Nose Bumper and RespiraGel Mustache, Baby Chin Strap, Baby Head Band and Circuit Bumpers, and Baby Bands for Chin and Head are trademarks of Global Respiratory Solutions, Inc. Babi.Plus is a registered trademark and GiO and Venti.Plus are trademarks of Galemmed Corporation. C1030.5