

Preemie Beenie™ Poly-Lined Knit Hat

respiralogics™

Poly-lined knit cap intended to decrease heat loss in pre-term and full-term infants

Preemie Beenie Poly-Lined Knit Hats are specifically designed to decrease heat loss in ELBW and Term infants in the critical early days of life. The Poly-Lined Knit Hat can significantly decrease cranial heat loss, where most evaporative heat loss occurs.

Our design provides a polyurethane liner between two layers of soft, knit fabric for comfort and prevention of heat loss. The Preemie Beenie knit maintains the shape of the Hat for a snug fit.


Ordering Information

45105	Preemie Beenie – Size 1 (≤ 25 weeks GSA) – White	50/box
45205	Preemie Beenie – Size 2 (25 to 28 weeks GSA) – Blue and white striped	50/box
45305	Preemie Beenie – Size 3 (≥ 29 weeks GSA) – Green and white striped	50/box

Companion Products

49050	StatStrap® Circuit Holder	50/box
-------	---------------------------	--------

Respiralogics products are available from Respiralogics and our specialty respiratory care distributors. For our latest product information visit www.respiralogics.com or email us at 4info@respiralogics.com.



FEATURES AND BENEFITS

Soft, stretchable knit for easy application

Snug fit to minimize movement and pressure on the head

Latex-free, DEHP-free polyurethane liner helps decrease heat loss in small infants

Color easily identifies Preemie Beenie size

Three sizes available to accommodate ELBW infants to full term

Single patient use, soft, skin-friendly two-ply knit hat

Respiralogics is a Global Respiratory Solutions, Inc. Company
3545 Airway Drive, Suite 104 • Reno, NV 89511 • 775 954 0160 • www.respiralogics.com

© 2024 Respiralogics. All rights reserved. Respiralogics, the Respiralogics logo, and Preemie Beenie are trademarks, and StatStrap is a registered trademark of Global Respiratory Solutions, Inc. C45000.2