

# StatStrap® Circuit Holder

Highly visible circuit holder clips to bedding, taking strain off tubes, circuits and cables

respiralogics™



The StatStrap Circuit Holder securely wraps around monitor cables and clips them to bed linens or patient gown to conveniently keep cables secure and tubing within easy reach.

The StatStrap Circuit Holder is designed for single patient use in neonatal critical care, labor and delivery, transport during Kangaroo Care and special procedures to safely secure tubing to clothing or bedding. The integral clip helps keep weight off cable connections for better monitor performance.

**Caution:** Not recommended for MRI use.

## Ordering Information

49005	StatStrap Circuit Holder	50/bag
49010	StatStrap Circuit Holder	100/box (2 bags of 50 ea)

Respiralogics products are available from Respiralogics and our specialty respiratory care distributors. For our latest product information visit [www.respiralogics.com](http://www.respiralogics.com) or email us at [4info@respiralogics.com](mailto:4info@respiralogics.com).



## FEATURES AND BENEFITS

Keeps circuit and cables stable with improved patient movement

Easy attachment to bed clothing without damaging tubing

Helpful during Kangaroo Care

Protects circuits, cables and tubes from damage

Strong Velcro® fastener and metal clip opens and closes without loss of integrity

Great for ambulatory patients

Single patient use, latex-free and DEHP-free

Respiralogics is a Global Respiratory Solutions, Inc. Company  
3545 Airway Drive, Suite 104 • Reno, NV 89511 • 775 954 0160 • [www.respiralogics.com](http://www.respiralogics.com)

© 2024 Respiralogics. All rights reserved. Respiralogics and the Respiralogics logo are trademarks, and StatStrap is a registered trademark of Global Respiratory Solutions, Inc. C49000.1