

# Jimmy Proof Straps™

Ventilator circuit security strap to prevent accidental disconnect

respiralogics™

Jimmy Proof Straps are an anti-disconnect device intended to secure the ventilator circuit connector to the tracheostomy tube holder to help reduce risk of accidental disconnects. The Straps are made of cotton with reinforced stitching to prevent stretch and rolling of the strap. The ventilator connector threads through the center opening of the Strap and the tapered ends of the Strap attach easily to the tracheostomy tube collar. Engineered to hold its shape with hook tab connectors to ensure quick release when disconnect is indicated.

The Jimmy Proof Straps are intended for use in critical care units, special procedure units as well as alternate care facilities and home care.

Jimmy Proof Straps were developed with special care by a dad wanting to make a difference for his son and change the quality of life for all tracheostomy patients. His engineering expertise and commitment to finding a better solution is found in the Jimmy Proof Straps and Danny Ties™.


## Ordering Information

DWL105005	Jimmy Proof Straps	250/case
DWL105005B	Jimmy Proof Straps	25/box

## Companion Products

DWL105001	Danny Ties, Small
DWL105002	Danny Ties, Medium
DWL105003	Danny Ties, Large
SP02	Sil.Flex Stoma Pad, Size 2
SP03	Sil.Flex Stoma Pad, Size 3
DP00	Sil.Flex TC Pad, 5 mm
DP01	Sil.Flex TC Pad, 10 mm

Shown with universal adaptor  
(universal adaptor not sold  
with Jimmy Proof Strap)



Respiralogics products are available from Respiralogics and our specialty respiratory care distributors. For our latest product information visit [www.respiralogics.com](http://www.respiralogics.com) or email us at [4info@respiralogics.com](mailto:4info@respiralogics.com).



## FEATURES AND BENEFITS

One piece design, constructed of soft cotton with reinforced stitched edges preventing stretching and rolling

Strengthened center hole of the Strap provides a stable opening for ventilator connector

Unique tapered, hook tabs easily adhere to Danny Ties and many commercially available tracheostomy collars

Cotton material absorbs moisture and wicks away from patient

Single patient use item, individually packed in environmentally friendly packaging

Respiralogics is a Global Respiratory Solutions, Inc. Company  
3545 Airway Drive, Suite 104 • Reno, NV 89511 • 775 954 0160 • [www.respiralogics.com](http://www.respiralogics.com)

© 2019 Respiralogics. All rights reserved. Respiralogics and the Respiralogics logo are trademarks of Global Respiratory Solutions, Inc. Sil.Flex is a trademark of Galemed Corporation. Danny Ties and Jimmy Proof Straps are trademarks of Leckie Medical Products, Inc. CDWL105005.4