

# Preemie Beenie™ Poly-Lined Knit Hat

respiralogics™

Poly-lined knit hat intended to decrease heat loss  
in pre-term and full-term infants

Preemie Beenie Poly-Lined Knit Hats are specifically designed to decrease heat loss in ELBW and Term infants in the critical early days of life. The Poly-Lined Knit Hat can significantly decrease cranial heat loss, where most evaporative heat loss occurs.

Our design provides a polyurethane liner between two layers of soft, knit fabric for comfort and prevention of heat loss. The Preemie Beenie knit maintains the shape of the Hat for a snug fit.


## Ordering Information

45105	Preemie Beenie – Size 1 (≤ 25 weeks GSA) – Blue Stripes	50/box
45205	Preemie Beenie – Size 2 (25 to 28 weeks GSA) – White	50/box
45305	Preemie Beenie – Size 3 (≥ 29 weeks GSA) – Pink and Blue Stripes	50/box

## Companion Products

49050	StatStrap® Circuit Holder	50/box
-------	---------------------------	--------

Respiralogics products are available from Respiralogics and our specialty respiratory care distributors. For our latest product information visit [www.respiralogics.com](http://www.respiralogics.com) or email us at [4info@respiralogics.com](mailto:4info@respiralogics.com).



## FEATURES AND BENEFITS

Soft, stretchable knit  
for easy application

Snug fit to minimize  
movement and pressure  
on the head

Latex-free, DEHP-free  
polyurethane liner  
helps decrease heat loss  
in small infants

Color easily identifies  
Preemie Beenie size

Three sizes available  
to accommodate ELBW  
infants to full term

Single patient use,  
soft, skin-friendly  
two-ply knit hat

Respiralogics is a Global Respiratory Solutions, Inc. Company  
3545 Airway Drive, Suite 104 • Reno, NV 89511 • 775 954 0160 • [www.respiralogics.com](http://www.respiralogics.com)

© 2025 Respiralogics. All rights reserved. Respiralogics, the Respiralogics logo, and Preemie Beenie are trademarks, and StatStrap is a registered trademark of Global Respiratory Solutions, Inc. C45000.4