

Baby Cap™ and Circuit Bumpers™

For the baby who just wants a soft, comfortable cap

respiralogics™

Baby Cap is intended to reduce heat loss and assists securing nasal interface during nCPAP or NIV therapy. Baby Cap and Circuit Bumpers permit the prongs and circuit to be maintained in the best position for a relaxed fit.

Baby Cap is a single patient use, soft knit Cap with an opening at the crown. Seven sizes are available to accommodate infants weighing ≤ 10 kg. Each size is color coded with fringe around the crown opening. Single patient use item, latex free.

Baby Chin Strap may be used in conjunction with Baby Cap. The Baby Chin Strap offers support with a soft, skin-friendly strap that cups the chin and is secured to the Baby Cap with hook and loop tabs to gently close the mouth.

The Baby Cap and Circuit Bumpers work in conjunction with neonatal/pediatric ventilator circuits and oxygen therapy tubing. Designed for use in neonatal critical care units, delivery rooms, pediatric critical care units and special procedure units for stable positioning of nasal interface.



Ordering Information

| | | |
|-------|---------------------------------------|----------|
| 30301 | Baby Cap and Circuit Bumpers – Size 1 | 12 / box |
| 30302 | Baby Cap and Circuit Bumpers – Size 2 | 12 / box |
| 30303 | Baby Cap and Circuit Bumpers – Size 3 | 12 / box |
| 30304 | Baby Cap and Circuit Bumpers – Size 4 | 12 / box |
| 30305 | Baby Cap and Circuit Bumpers – Size 5 | 12 / box |
| 30306 | Baby Cap and Circuit Bumpers – Size 6 | 6 / box |
| 30307 | Baby Cap and Circuit Bumpers – Size 7 | 6 / box |

Companion Products


| | | | |
|-------|--|------|-------------------------------------|
| 1002 | Babi.Plus® nCPAP Kit | B020 | Babi.Plus nCPAP Nasal Cannula, Sz 0 |
| 1010 | Babi.Plus nCPAP Kit, Small | B021 | Babi.Plus nCPAP Nasal Cannula, Sz 1 |
| 1020 | Babi.Plus nCPAP Kit, Large | B022 | Babi.Plus nCPAP Nasal Cannula, Sz 2 |
| 1005 | Venti.Plus™ Single Pole Clamp | B023 | Babi.Plus nCPAP Nasal Cannula, Sz 3 |
| 1006 | Venti.Plus Dual Pole Clamp | B024 | Babi.Plus nCPAP Nasal Cannula, Sz 4 |
| 30100 | Baby Nose Bumper™ and RespiraGel Mustache™ | B025 | Babi.Plus nCPAP Nasal Cannula, Sz 5 |
| 30400 | Baby Chin Strap™ | B026 | Babi.Plus nCPAP Nasal Cannula, Sz 6 |
| | | B027 | Babi.Plus nCPAP Nasal Cannula, Sz 7 |

Respiralogics products are available from Respiralogics and our specialty respiratory care distributors. For our latest product information visit www.respiralogics.com or email us at 4info@respiralogics.com.



Respiralogics is a Global Respiratory Solutions, Inc. Company
3545 Airway Drive, Suite 104 • Reno, NV 89511 • 775 954 0160 • www.respiralogics.com

© 2019 Respiralogics. All rights reserved. Respiralogics, the Respiralogics logo, Baby Cap and Circuit Bumpers, Baby Nose Bumper and RespiraGel Mustache, Baby Head Band and Circuit Bumpers and Baby Chin Strap and trademarks of Global Respiratory Solutions, Inc. Venti.Plus is a trademark and Babi.Plus is a registered trademark of GaleMed Inc. C30300.1



FEATURES AND BENEFITS

Fashioned to provide a comfortable fit for all infants

Stretchable knit allows easy application

Circuit Bumpers attach firmly to the Cap and hold the inspiratory and expiratory limbs of the circuit in place

Hook strips are sewn into the Cap to prevent movement of the Circuit Bumpers

Allows caregivers to freely move infant without interrupting CPAP therapy

Crown opening for checking fontanels and accommodating scalp vein IVs

Useful in regulating temperature for infants

Single patient use, latex free