

# Baby Nose Bumper™ and RespiraGel Mustache™

For babies who want a pretty nose after therapy

respiralogics™

Baby Nose Bumper and RespiraGel Mustache provide optimal fixation of nasal prongs for infants during nCPAP or NIV therapy. The mustache is made of RespiraGel™, our new “skin friendly” hydrocolloid-based adhesive. Baby Nose Bumper wraps around the body of the nasal prong and gently holds the nasal interface to the mustache, providing a secure grip and gentle cushion for the nares. Baby Nose Bumper and RespiraGel Mustache are single-patient-use items.

Baby Nose Bumper and RespiraGel Mustache may be used in conjunction with Babi.Plus® Nasal Prongs and other commercially available non-invasive interfaces in neonatal critical care units, delivery rooms, pediatric critical care units and special procedure units.



## Ordering Information

30100	Baby Nose Bumper and RespiraGel Mustache	20 / box
30150	Baby Nose Bumper Replacement (25 bags of 4 each)	100 / box

## Companion Products

1002	Babi.Plus nCPAP Kit
1010	Babi.Plus nCPAP Kit, Small
1020	Babi.Plus nCPAP Kit, Large
30400	Baby Chin Strap™
B020 series	Babi.Plus nCPAP Nasal Cannula, Sizes 0 to 7
30500 series	Baby Head Band and Circuit Bumpers™, Sizes fitting 18 to 42 cm
30300 series	Baby Cap™ and Circuit Bumpers™, Sizes 1 to 7
30510	Baby Band™ for Chin and Head

Respiralogics products are available from Respiralogics and our specialty respiratory care distributors. For our latest product information visit [www.respiralogics.com](http://www.respiralogics.com) or email us at [4info@respiralogics.com](mailto:4info@respiralogics.com).



## FEATURES AND BENEFITS

Baby Nose Bumper and RespiraGel mustache provide a secure, stable fit of the prongs in the nares

The mustache is constructed of RespiraGel, a latex-free, skin-protective hydrocolloid adhesive

Baby Nose Bumper gently grasps the soft hook on the mustache to provide a secure, stable grip

Baby Nose Bumper fits Babi.Plus Nasal Prongs

Baby Nose Bumper and RespiraGel Mustache allow caregivers to move infant without interrupting therapy

Respiralogics is a Global Respiratory Solutions, Inc. Company  
3545 Airway Drive, Suite 104 • Reno, NV 89511 • 775 954 0160 • [www.respiralogics.com](http://www.respiralogics.com)

© 2019 Respiralogics. All rights reserved. Respiralogics, the Respiralogics logo, Baby Nose Bumper and RespiraGel Mustache, Baby Head Band and Circuit Bumpers, Baby Chin Strap, Baby Cap and Circuit Bumpers, and Baby Band for Chin and Head are trademarks of Global Respiratory Solutions, Inc. Babi.Plus is a registered trademark of GaleMed Inc. C30100.1