

Ories™ Cohesive Bandage

A soft wrap to secure catheters around tender skin

respiralogics™

Ories Cohesive Bandage is a non-adhesive wrap constructed from 100% polyester made by a special weaving method. The wrap adheres to itself, not to skin, hair or clothing. Ories material allows it to conform to the body while maintaining the position of a pad or catheter. Once the site is wrapped with Ories, there is no need for another tape to secure the wrap.



Ories is suitable for fixation of primary and secondary dressings, pads, padding material, catheters and cannulas. The material is easy to tear, breathable, stretchable, water repellant and soft to the touch for tender skin. The high stretchability makes for easy wrapping. The layers of Ories do not slip so it can securely attach to itself, not to skin.

Ories is a skin-friendly, hypoallergenic, single-patient-use wrap available in two sizes: 2.5 cm width and 3.75 cm width. Each cohesive bandage is individually packaged in clear plastic.

Ordering Information

FMMJ3001	Ories 2.5 cm width by 3.5 m length	25/box
FMMJ3002	Ories 3.75 cm width by 3.5 m length	18/box

Respiralogics products are available from Respiralogics and our specialty respiratory care distributors. For our latest product information visit www.respiralogics.com or email us at 4info@respiralogics.com.



FEATURES AND BENEFITS

- Tears easily by hand, no scissors needed
- Anti-slip material to prevent movement
- Non-adhesive material that securely attaches to itself
- Water repellant
- Breathable material to prevent skin irritation
- Soft-touch wrap prevents skin irritation
- High stretchability to conform to patient/application
- Single patient use

Respiralogics is a Global Respiratory Solutions, Inc. Company
3545 Airway Drive, Suite 104 • Reno, NV 89511 • 775 954 0160 • www.respiralogics.com

© 2026 Respiralogics. All rights reserved. Respiralogics and the Respiralogics logo are trademarks of Global Respiratory Solutions, Inc. Ories is a trademark of FOXXMED Ltd. FMMJ3000.3