

Venti.Plus™ Test Lung

Simulates respiratory system

respiralogics™

The Venti.Plus Test Lung is a great tool to simulate the respiratory system, demonstrate mechanical ventilator applications and perform ventilator circuit testing prior to clinical use.

The Venti.Plus Test Lung is compact in design, lightweight and portable. The Venti.Plus Test Lung – 1.0 Liter is packaged with a 1.0 Liter Latex Free Silicone Bag to simulate the respiratory system of an adult lung. The Venti.Plus Test Lung – 0.5 Liter is packaged with a 0.5 Liter Latex Free Silicone Bag to simulate the respiratory system of a pediatric lung. The convenient turn and lock design at the base of the lung allows for cleaning and decontamination of the Silicone Bag. Each Test Lung is calibrated for accuracy and validated for resistance and compliance in the application range to ensure compliance with specification.

The Venti.Plus Test Lung is the ideal instrument for training and demonstration in addition to performing ventilator verification testing.

Test Lung Specifications

	Adult	Child
Part Number	2000	2005
Size	1.0 Liter	0.5 Liter
Compliance	0.02 L/cm H ₂ O	0.01 L/cm H ₂ O
Resistance	Rp5 (reference)	Rp5 (reference)
Volume	0 – 600 ml	0 – 400 ml
Weight	230 grams	140 grams

Ordering Information

2000	Venti.Plus Test Lung – 1.0 L	1 each
2005	Venti.Plus Test Lung – 0.5 L	1 each

Respiralogics products are available from Respiralogics and our specialty respiratory care distributors. For our latest product information visit www.respiralogics.com or email us at 4info@respiralogics.com.



FEATURES AND BENEFITS

- Calibrated for higher degree of accuracy
- ISO compliant 15 mm connection port
- Lightweight, easy to use and portable
- Manufactured of high performance material
- Convenient turn and lock design for decontamination
- Packaged in convenient, reusable container

Respiralogics is a Global Respiratory Solutions, Inc. Company
3545 Airway Drive, Suite 104 • Reno, NV 89511 • 775 954 0160 • www.respiralogics.com

© 2019 Respiralogics. All rights reserved. Respiralogics and the Respiralogics logo are trademarks of Global Respiratory Solutions, Inc. Venti.Plus is a trademark of Galemed Corporation. CTL2000.5