

Babi.Plus® nMask

Silicone resuscitation mask for optimal comfort and seal

respiralogics™

Babi.Plus nMask is precisely designed to comfortably conform to an infant's face, facilitating a good seal for optimal ventilation and non-invasive support.



The circular design of the nMask allows easy creation of a comfortable anatomical seal for the very small premature infant to toddler face. The clear nMask provides a good view of the patient/mask interface to ensure proper fit and size.

The silicone nMask is available in six sizes; ELBW through small pediatric patient. Each nMask is individually packaged and ready for use with resuscitation bag or T-piece resuscitator system.

The Babi.Plus nMask can be used in all neonatal critical care units, delivery rooms, pediatric critical care units and special procedure units.

Ordering Information

AN0038	nMask 25 mm	1 piece/bag, 5 bags/box
AN0039	nMask 30 mm	1 piece/bag, 5 bags/box
AN0040	nMask 35 mm	1 piece/bag, 5 bags/box
AN0041	nMask 40 mm	1 piece/bag, 5 bags/box
AN0042	nMask 45 mm	1 piece/bag, 5 bags/box
AN0043	nMask 50 mm	1 piece/bag, 5 bags/box

Respiralogics products are available from Respiralogics and our specialty respiratory care distributors. For our latest product information visit www.respiralogics.com or email us at 4info@respiralogics.com.



FEATURES AND BENEFITS

- Designed specifically to conform to premature and newborn faces
- Ergonomic design with soft material to provide better seal and comfort
- Silicone material to minimize skin irritation; bio-compatible, hypoallergenic and latex free
- Available in six sizes; fits ELBW infant through small pediatric faces
- 15mm connector with two-level barb design for easier attachment to resuscitator devices or T-piece resuscitator circuit
- Single patient use
- Individually packaged and easy to apply to neck of resuscitation instrument

Respiralogics is a Global Respiratory Solutions, Inc. Company
3545 Airway Drive, Suite 104 • Reno, NV 89511 • 775 954 0160 • www.respiralogics.com

© 2019 Respiralogics. All rights reserved. Respiralogics and the Respiralogics logo are trademarks of Global Respiratory Solutions, Inc. Babi.Plus is a registered trademark of Galemed Corporation. CAN0038.2