


For Immediate Release

**Respiralogsics introduces 30500 Series Baby
Head Band™ and Circuit Bumpers™**
Reno, NV - Spring 2018



Introducing 30500 Series Baby Head Band™ and Circuit Bumpers™

Finally, a solution for clinicians needing a more comfortable and versatile way to secure nasal prongs on their most demanding patients.

Baby Head Band and Circuit Bumpers are intended to secure the nasal interface during oxygen therapy, nCPAP or NIV therapy.

The Baby Head Band is made of breathable, non-latex foam that absorbs and wicks away moisture leaving the scalp and skin dry. The soft material of the head band allows it to conform to even the tiniest of heads while maintaining the circuit position. The soft outer covering of the Baby Head Band allows for easy head movement while holding the nasal prongs and circuit stable.

The new reversible Circuit Bumpers secure and stabilize the position of the inspiratory and expiratory limbs to the soft outer side of Head Band while providing optimal positioning for infants of all shapes and sizes.

The Head Band and Circuit Bumpers work in conjunction with neonatal/pediatric ventilator circuits and oxygen therapy tubing.

Benefits:

- Fashioned to provide a comfortable fit for all infants
- Sized by head circumference for easy and quick application
- Reversible Velcro® Bumpers allow for secure attachment and minimize pressure points
- Allows caregivers to freely move infant without interrupting therapy
- Opened at the crown for fontanel checks and to accommodate scalp vein IV's
- Available in four (4) sizes for very low birth weight up to 10 kg infants
- Single Patient Use, non-latex material

Designed for use in NICU, PICU, surgery and clinics for stable positioning of nasal interface.

Respiralogsics specializes in providing innovative respiratory care products to meet the challenging needs of respiratory patients. To benefit clinicians, we have joined with leading specialty respiratory distributors to provide expert local assistance and the training required for successful implementation of our products.

The Respiralogsics complete line of products are available from Respiralogsics and our specialty respiratory care distributors. Visit www.respiralogsics.com for all Respiralogsics products and information.

Respiralogsics, Baby Head Band and Circuit Bumpers are trademarks of Global Respiratory Solutions, Inc.

For more information, contact Respirologics at 1.775.954.0160

Company: Respirologics

Address: 3545 Airway Drive, Suite 104, Reno, NV 89511

Telephone: 1.775.954.0160

Web: www.respirologics.com

Contact: Beth Keifer

respirologics[™]

A Global Respiratory Solutions, Inc. Company